

Shamir Autograph III SV™


A highly personalized single vision lens

What is it

Autograph III SV™ is a premium freeform single vision lens designed with EyePoint Technology III to incorporate power corrections based on the precise measurement of customers' individual prescriptions, and simulation of their perceived viewing fields. Patients wearing Autograph III SV™ will experience the sharpest and most comfortable vision available. It does this by using the highest degree of personalization and by precisely matching the viewing experience with any frame choice.

As-worn Quadro technology is also part of the design, to guarantee that visual comfort and acuity give your patient the optimal vision possible.


Autograph III SV™ Lens Markings



Benefits

Optimal visual experience
with any frame of choice

The widest possible
fields of clear vision



Autograph III SV: Fitting Height is required Available with Glacier Plus or Glacier Anti-Fog premium coatings*	Range [D]	Cylinder
1.5: Clear, Polarized, Transitions Gen 8, XTRActive™	-12.00 to +6.00	to -6.00
1.53 Trivex: Clear, Blue Zero, Polarized, Transitions Gen 8	-14.00 to +6.75	to -6.50
1.59 Poly: Clear, Blue Zero, Polarized, Transitions Gen 8, XTRActive™	-14.00 to +8.75	to -7.00
1.6: Clear, Blue Zero, Polarized, Transitions Gen 8, XTRActive™	-14.00 to +8.00	to -7.00
1.67: Clear, Blue Zero, Polarized, Transitions Gen 8, XTRActive™	-16.00 to +16.00	to -8.00
1.74: Clear, Transitions Grey	-20.00 to +15.50	to -9.00

*Glacier Anti-Fog coating is available on Trivex, Poly and 1.6

Featured Technologies



Shamir Canada 1828 Blue Heron Drive, Unit #37, London, ON N6H 0B7
TF 1-844-710-5367 | T 1-519-471-0750 | F 1.855.742.6471 | www.shamirlens.ca